

**Product number: R-PE**

**Product name: Lyophilized R-Phycoerythrin**

## General Data

**Molecular Mass:** 240 kDa

**Solubility:** Water, Aqueous Buffers

**Insoluble:** Acetone, Chloroform, Toluene

**Storage:** Store in absence of light, desiccate and refrigerate. Do not freeze.

## Description

Lyophilized **R-PE** is a lyophilized phycobiliprotein from water with sugar as additive. No ammonium sulfate or other materials that may interfere with your conjugation process are added to this product. **R-PE** consists of  $\alpha$ ,  $\beta$  and  $\gamma$  subunits and is present as  $(\alpha\beta)_6\gamma$ .

## Applications

- Immunoblotting
- Immunostaining
- Resonance Energy Transfer (RET) Acceptor
- Flow Cytometry

## Advantages


- Perfectly suited for excitation with the 488-nm or 532-nm diode lasers
- Sensitive; high extinction coefficient and high quantum yield
- Good aqueous solubility; this label does not alter the solubility of bioconjugates

## Spectral Data

**Solvent System:** phosphate buffer pH 7.4

Sample	$A_{565}/A_{280}$	$A_{565}/A_{498}$	$A_{620}/A_{565}$	Absorption max. [nm]	Extinction Coefficient [ $M^{-1}cm^{-1}$ ]	Emission <sup>1</sup> max. [nm]
R-PE	> 5.4	> 1.5	> 0.005	565, 539, 498	1,960,000 at 565 nm	573

<sup>1</sup> Excitation at 525 nm



Absorption and emission spectrum of **R-PE** in phosphate buffer (pH 7.4)

**Product number: R-PE**

**Product name: Lyophilized R-Phycoerythrin**

---

## Important Notes

**Reconstitution:** Reconstitute the whole bottle of lyophilized **R-PE** (5 mg) with your specific volume of conjugation buffer to adjust the concentration for further use.

**Weight:** One bottle of lyophilized **R-PE** contains about 5 mg of **R-PE** with added sugar as protective. Please refrain from determining the **R-PE** concentration directly by weight. In order to obtain an accurate concentration of **R-PE** in mg/mL, use the extinction coefficient and determine the concentration using the following formula:

$$[\mathbf{R-PE}] = 0.122 \times A_{565}$$

where **[R-PE]** is the concentration of **R-PE** in mg/mL and  $A_{565}$  is the absorbance at 565 nm, provided  $A_{565}$  is in the range between 0.3 to 0.8

**Usage:** No preservative ( $\text{NaN}_3$ ) is added to the product. Once the protein is reconstituted it should be used as soon as possible.